



Walletz4u®

Nitrile Exam Gloves Powder Free, Standard Cuff

Walletz4u®, the industry's lowest-weight Accelerator-Free Premium Nitrile Examination Gloves that is compliant with both US ASTM and European EN regulatory requirements. As such, Walletz4u® will be able to offer a "Low Dermatitis Potential" labeling claim, as supported by clinical testing. Given the unique formulation used to manufacture it, Walletz4u® is also sulphur-free.



Walletz4u®		
Length (mm)	≥ 230	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.07 ± 0.02	
Finger (13mm ± 3mm from tip)	0.09 ± 0.02	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 400
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.5
Physical Dimensions	S2	4.0
Physical Properties	S2	4.0
Visual Inspection (Major)	S4	2.5
Visual Inspection (Minor)	S4	4.0
Particulate Residue	N = 5	≤ 2mg/glove

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Sapphire Blue colour

PACKAGING

100 gloves per box
10 boxes per carton

REGULATORY COMPLIANCE

FDA 510(k), MDD 93/42/EEC,

STANDARDS

ASTM D6319, ASTM D412, ASTM D573,
ASTM D5151, ASTM D6124,
EN455 part 1, 2, & 3,

MANUFACTURING ACCREDITATIONS

ISO 9001:2015
ISO 13485:2016
ISO 14001:2015
ISO 45001:2018

Chemotherapy Drugs and Concentration
(Tested for Resistance to Permeation by Chemotherapy
Drugs as per ASTM D6978-05)

Minimum Breakthrough
Detection Time (minutes)

Carmustine (BCNU), 3.5mg/ml (5,300 ppm)	Not recommended
Cisplatin, 1.0mg/ml (1,000 ppm)	>240 minutes
Cyclophosphamide (Cytosan), 20.0mg/ml (20,000 ppm)	>240 minutes
Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm)	>240 minutes
Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm)	>240 minutes
Etoposide (Toposar), 20.00mg/ml (20,000 ppm)	>240 minutes
Fluorouracil, 50.0mg/ml (50,000 ppm)	>240 minutes
Methotrexate, 25.0mg/ml (25,000 ppm)	>240 minutes
Mitomycin C, 0.5mg/ml (500 ppm)	>240 minutes
Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm)	>240 minutes
Thiotepa, 10.0mg/ml (10,000 ppm)	Not recommended
Vincristine Sulfate, 1.0mg/ml (1,000 ppm)	>240 minutes

WARNING: Gloves used for protection against chemotherapy drug exposure should be selected specifically for the type of chemicals being used. Due to the variety and concentration of chemotherapy drugs used in treatments, the resistance table shown does neither warrant nor imply the safe use of the gloves against chemotherapy drugs resistance in every case. The safe use of gloves in chemotherapy treatment is solely the decision of clinicians authorised to make such decision.



Reg. No: GMD68412486918A

Walletz4u
Manufactured by:
J&B ENGINEERING & MEDICAL TECHNOLOGIES (M) SDN BHD

A brand by
Saw & Meam

Marketed & Distributed by
Goldmark

MADE IN MALAYSIA